

**IB DP Chemistry
Prescribed practicals
lab manual**

www.msjchem.com

Prescribed practicals

Included in the “Applications and skills” sections of the guide are a series of practicals that students must cover either in a laboratory environment or as a simulation. The skills and general techniques associated with these common practicals may be assessed as part of the external assessment. The list of required practicals for chemistry is shown in the table below with suggested practicals.

Topic 1.2	Obtaining and using experimental data for deriving empirical formulas from reactions involving mass changes.
	Determine the formula of MgO Determine the formula of hydrated copper sulfate
Topic 1.3	Use of the experimental method of titration to calculate the concentration of a solution by reference to a standard solution.
	Preparation of a standard solution/acid base titration Determine the percentage CaCO ₃ in eggshells
Topic 1.3	Obtaining and using experimental values to calculate the molar mass of a gas from the ideal gas equation.
	Determine the molar mass of a gas (CO ₂) experimentally
Topic 5.1	A calorimetry experiment for an enthalpy of reaction should be covered and the results evaluated.
	Calculate enthalpy change of neutralisation (HCl and NaOH) Calculate enthalpy change Zn and CuSO ₄
Topic 6.1	Investigation of rates of reaction experimentally and evaluation of results.
	Sodium thiosulfate and hydrochloric acid reaction Iodine clock reaction
Topic 8.2	Candidates should have experience of acid–base titrations with different indicators.
	Covered in topic 1.3
Topic 8.3	Students should be familiar with the use of a pH meter and universal indicator.
	Reactions of acid and bases lab
Topic 9.2	Performance of laboratory experiments involving a typical voltaic cell using two metal/metal–ion half-cells.
	Voltaic cell Zn and Cu half cells
Topic 10.1	Construction of 3D models (real or virtual) of organic molecules.
	Molecular shapes simulation https://phet.colorado.edu/en/simulation/molecule-shapes Molymods http://www.molymod.com/ Chemdoodle https://www.chemdoodle.com/
Topic 19.1	Perform lab experiments which could include single replacement reactions in aqueous solutions.
	Displacement reactions

Topic 1.2 - Obtaining and using experimental data for deriving empirical formulas from reactions involving mass changes.

Title: Determine the empirical formula of MgO experimentally


Background: When magnesium is heated in air, it reacts with oxygen. During this oxidation reaction, magnesium oxide is produced. This increases the mass. If we know the mass of magnesium at the start, and the mass of magnesium oxide produced at the end, we can calculate the mass of oxygen which has been combined with the magnesium. We can use these masses to determine the empirical formula of magnesium oxide.

Procedure:

1. Record the mass of the empty crucible with its lid and write down the result in a table.
2. Clean the piece of magnesium ribbon with sandpaper, then coil it loosely around a pencil.
3. Put the magnesium ribbon into the crucible and put the lid on. Record the mass of the crucible, lid, and magnesium together, and write down the result in your table.
4. Put the crucible onto the pipe-clay triangle. Leave the lid slightly ajar to allow air into the crucible. Heat gently for a minute, then strongly. Continue heating until the reaction has finished (the magnesium will glow at first, then look a bit like a furry grey-black caterpillar)
5. Turn the gas burner off and allow the crucible to cool for a few minutes. Record the mass of the crucible with its lid and contents.

Results:

	Mass ($\pm 0.01\text{g}$)
crucible + lid	
crucible + lid + magnesium	
crucible + lid + contents after reaction	



Calculations


1. Calculate the mass of magnesium used in the reaction.
2. Calculate the mass of oxygen that combined with the magnesium.
3. Convert the mass of magnesium and oxygen to amount (in mol).
4. Find the lowest whole number ratio of Mg to O – this will give the empirical formula.

Title: Determine the formula of hydrated copper(II) sulfate

Background: Many crystals occur commonly in hydrated form, which means that they contain water molecules within the crystal structure in definite proportions. This water is called water of crystallization. The formula for these crystals shows the number of water molecules present per formula unit of crystal using a dot before the water. When these crystals are heated, they lose their water of crystallization and are then said to be anhydrous.

Procedure:

1. Weigh the empty crucible, and then weigh into it between 2.00 g and 3.00 g of hydrated copper(II) sulfate. Record all masses accurate to the nearest 0.01 g.
2. Support the crucible securely in the pipe-clay triangle on the tripod over the Bunsen burner.



3. Heat the crucible and contents, gently at first, over a medium Bunsen flame, so that the water of crystallisation is driven off steadily. The blue colour of the hydrated compound should gradually fade to the greyish-white of anhydrous copper(II) sulfate. **Avoid over-heating, which may cause further decomposition, and stop heating immediately if the colour starts to blacken. If over-heated, toxic or corrosive fumes may be evolved.**
4. Allow the crucible and contents to cool. Re-weigh the crucible and contents once cold.
5. 'Heat to constant mass' by repeating steps (iii) and (iv) until you get consistent readings.

Results:

	Mass (± 0.01 g)	Observations
Mass of crucible + lid		
Mass of crucible + lid + $\text{CuSO}_4 \cdot x\text{H}_2\text{O}$ before heating		
Mass of crucible + lid + CuSO_4 (anhydrous) after heating		

Calculations:

1. Calculate the molar masses of H_2O and CuSO_4 (Relative atomic masses: H=1.01, O=16.00, S=32.06, Cu=63.55).
2. Calculate the mass of water driven off, and the mass of anhydrous copper(II) sulfate formed in your experiment.
3. Calculate the amount (in mol) of anhydrous copper(II) sulfate formed
4. Calculate the amount (in mol) of water driven off.
5. Calculate amount (in mol) of water would have been driven off if one mole of anhydrous copper(II) sulfate had been formed.
6. Determine the formula for hydrated copper(II) sulfate.

Topic 1.3 - Use of the experimental method of titration to calculate the concentration of a solution by reference to a standard solution.

Title: Preparation of a standard solution

Background: Standard solutions are solutions with known concentrations, generally to four significant figures. There are two different ways to make a standard solution. We can make a primary or a secondary standard. A primary standard is prepared directly by dissolving a known mass of sample to make a known volume of solution. A secondary standard is prepared by dissolving an approximate amount of sample into a volume of solvent and determining its exact concentration through titration experiments. Primary standards are prepared from compounds that are at least **99.9% pure**, have a **definite composition**, are **water soluble**, are **easily weighed**, and **do not change composition on contact with air**. Oxalic acid dihydrate ($\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$) fits these criteria and therefore may be used as a primary standard. Sodium hydroxide absorbs water when it comes into contact with air and therefore it is difficult to obtain a pure, dry sample to weigh. For this reason, the sodium hydroxide solution will be titrated with the oxalic acid standard to become a secondary standard.

In the first part of this experiment you will prepare a solution of known concentration of oxalic acid. The oxalic acid crystallizes with two water molecules per oxalic acid in the crystalline network. For this reason, we will weigh out an appropriate amount of oxalic acid dehydrate to dissolve in water. The water molecules in the crystal network will become part of the water of solution once it is dissolved. The molar concentration of oxalic acid dihydrate will be the same as the molar concentration of oxalic acid.

Preparation of primary standard

Procedure:

1. Obtain a 100 cm³ beaker.
2. Measure between 1.25 g and 1.30 g of pure oxalic acid dehydrate crystals and place into the beaker.
3. Add approximately 30 cm³ of deionized water to the beaker and dissolve the crystals.
4. Transfer the solution into a clean 100.00 cm³ volumetric flask.
5. Rinse the beaker with 20 cm³ of deionized water and pour this solution into the volumetric flask and repeat. This will ensure that all of the oxalic acid is transferred into the volumetric flask.
6. Fill the volumetric flask to within about 2 cm³ of the mark and allow it to sit for a minute. This will allow any water clinging to the edges of the neck to drain into the flask.
7. Using a pipette, fill the flask to the 100.00 cm³ mark with water.
8. Stopper the flask and mix the solution by repeated inversion and swirling.


Calculations:

Mass of oxalic acid =

1. Calculate amount (in mol) of oxalic acid
2. Calculate the concentration of the oxalic acid solution using the equation $C = n/V$.

Standardisation of sodium hydroxide solution

In order to standardize the sodium hydroxide solution, you will perform a titration. Sodium hydroxide reacts with oxalic acid according to the reaction below:



You will measure a 20.00 cm³ aliquot (an aliquot is an exact quantity of a substance or solution) of the oxalic acid solution into a flask and add an indicator. An indicator is a substance that changes colour when a solution changes from acidic to basic. The common indicator used for acid base titrations is phenolphthalein. Phenolphthalein is colourless in a solution that is acidic and bright pink in a solution that is basic. In this titration the oxalic acid solution is acidic and therefore phenolphthalein will be colourless. The sodium hydroxide solution will be added drop wise from a burette into the flask containing the oxalic acid and indicator. As the sodium hydroxide is added to the flask it will react with the oxalic acid and be neutralized. At the point where all of the oxalic acid is reacted, the next drop of sodium hydroxide will make the entire solution basic and it will turn pink. At this point you have completed the titration. In order to get the best precision possible, you should repeat each titration until you get 3 trials that are within 0.10 cm³ of each other (concordant).

Results

	Trial 1 (rough)	Trial 2	Trial 3	Trial 4
Volume of oxalic acid solution ± 0.08 cm ³	20.00	20.00	20.00	20.00
NaOH buret, initial volume reading ± 0.05 cm ³				
NaOH buret, final reading ± 0.05 cm ³				
Volume of NaOH solution used ± 0.10 cm ³				

Calculations:

1. Calculate the average volume of NaOH solution used (choose values within 0.10 cm³)
2. Calculate the amount (in mol) of oxalic acid in 20 cm³ of oxalic acid solution.
3. Use the ratio of NaOH to H₂C₂O₄ to find the amount (in mol) of NaOH.
4. Use the equation $C=n/V$ to find concentration of NaOH solution

Title: Determine the percentage of calcium carbonate in egg shells

Background: In this experiment, you will determine the percentage of calcium carbonate in an egg shell. Calcium carbonate reacts with hydrochloric acid according to the following equation:



It is difficult to determine exactly when all the solid calcium carbonate has reacted to one drop of standard hydrochloric acid solution if it is titrated directly. In this experiment a known amount of excess acid is added to the sample to ensure that all the calcium carbonate has reacted. The excess acid is then diluted and made up to a known volume. Aliquots of this diluted excess acid solution are then titrated with a standard solution of sodium hydroxide.

Procedure:

1. Clean and dry an eggshell, removing as much of the egg membrane as possible.
2. Weigh 0.60 g of crushed eggshell and place in a conical flask.
3. Slowly add 20.00 cm³ of 1.00 mol dm⁻³ HCl to the egg shell in the conical flask.
4. When the reaction has finished, add about 20 cm³ of distilled water to the solution.
5. Transfer the solution to a 100.00 cm³ volumetric flask and make up to 100.00 cm³ with distilled water.
6. Measure 10.00 cm³ of this solution and place in a clean conical flask.
7. Add a few drops of phenolphthalein indicator and titrate with 0.100 mol dm⁻³ NaOH.

Results:

Trial	1 (rough)	2	3	4	5
Initial volume of NaOH (±0.05 cm ³)					
Final volume NaOH (±0.05 cm ³)					
Total volume NaOH used (±0.10 cm ³)					
Colour change					

Choose results that are within 0.10 cm³ of each other and calculate the average volume of NaOH used.

Average volume of NaOH =

Calculations:

1. Calculate the amount (in mol) of hydrochloric acid present initially in 20.00 cm³ of 1.00 mol dm⁻³ HCl_(aq)
2. Calculate the amount of sodium hydroxide (in mol) present in the volume of 0.100 mol dm⁻³ NaOH_(aq) used for the titration.
3. Determine the amount (in mol) of unreacted HCl_(aq) present in the 10.00 cm³ sample.
4. Determine the amount (in mol) of unreacted HCl_(aq) present in the 100 cm³ solution.
5. Determine the amount (in mol) of HCl_(aq) that reacted with the egg shell.
6. Calculate the amount (in mol) of calcium carbonate present in the egg shell.
7. Calculate the mass of calcium carbonate present in the egg shell.
8. Determine the percentage by mass of calcium carbonate in the egg shell.
9. Determine the total percentage uncertainty associated with your result.

Topic 1.3 Obtaining and using experimental values to calculate the molar mass of a gas from the ideal gas equation.

Title: Determine the molar mass of CO₂ experimentally

Background: In this lab, we are going to determine the molar mass of CO₂ experimentally. A simple calculation using a periodic table would give us the correct answer for the molar mass of CO₂, however, you are going to conduct an experiment in the lab to see how close you can come to the accepted value. The ideal gas law equation ($PV = nRT$) can be rearranged to solve for n , the amount in mol:

$$n = \frac{PV}{RT}$$

Molar mass (M), is expressed in units of grams per mole (g mol^{-1}). The mass of the CO₂ released from the pressurized container the volume of the gas will be measured. You can use the method of water displacement at room conditions and then substitute the measurements of volume, temperature, and pressure into the ideal gas law equation in order to find n , the amount (in mol) of CO₂. Once we know the mass, and the amount (in mol), we will divide the mass by the number of moles to calculate the molar mass (g mol^{-1}).

Procedure:

1. Determine the initial mass for the CO₂ canister provided by your instructor ($\pm .01$ gram).
2. Set up the water basin for collecting a gas in an inverted graduated cylinder by water displacement. Use a 100 cm³ graduated cylinder.
3. Release CO₂ from the pressurized container and collect a sample of gas with a volume of approximately 70.0 cm³ to 90.0 cm³. Collect every bubble that leaves the canister.
4. Adjust the cylinder up or down so that the water level inside the graduated cylinder is the same as the level outside the cylinder. **THIS IS VERY IMPORTANT!** If you miss this simple step, the pressures inside and outside the cylinder will not be equal, and you will have an inaccurate estimation of the pressure inside the cylinder where the gas is collected.
5. Record the volume of CO₂ that was collected in the Observations and Data section.
6. Measure and record the temperature of the water in the basin (± 0.1 °C). This will be used to determine the vapour pressure of water H₂O_(g) in the graduated cylinder.
7. Read the air pressure in the lab and convert this to Pa.
8. Measure and record the mass of the CO₂ canister after your experiment. (± 0.01 g).

Results and calculations:

Mass

Initial mass of CO₂ container (m_{initial}) = _____

Final mass of CO₂ container (m_{final}) = _____

Mass of CO₂ collected ($m_{\text{initial}} - m_{\text{final}}$) = _____

Volume

Volume of the gas collected (cm³) = _____

Volume of CO₂, converted to m³ (cm³ ÷ 10⁶) = _____

Temperature

Temperature of water in basin = _____

Temperature of water in kelvin (K) = _____

(assume that the water and CO₂ are the same temperature)

Pressure

Air pressure in the room = _____

Vapour pressure of water = _____

$P_{\text{total}} = P_{\text{CO}_2} + P_{\text{water vapour}}$ _____


Pressure of dry CO₂ _____

Calculate the molar mass of the gas using the equation:

$$M = \frac{mRT}{PV}$$

Title: Determine the enthalpy change for reaction of Zn and CuSO₄

Background: Chemical reactions can be classified as exothermic or endothermic depending on whether they have a negative or positive enthalpy change. In this experiment, you will calculate the enthalpy change for the following reaction:



Procedure:

1. Using a 50.0 cm³ measuring cylinder, measure 50.0 cm³ of 0.500 mol dm⁻³ CuSO₄ and pour into an insulated cup.
2. Using a mass balance, measure 3.00g of Zn powder.
3. Insert a temperature probe into the CuSO₄, start the data collection, and record the temperature every 30 seconds.
4. At 3 minutes (180 seconds), add the 3.00 g of Zn powder to the insulated cup containing the CuSO₄, and stir gently.
5. Continue to record the temperature every 30 seconds for a further 6 minutes (360 seconds).


Results:

Time (min)	0.0	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5
Temperature (°C)							not done			
Time (min)	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5
Temperature (°C)										

Analysis:

Plot a graph of temperature (y axis) against time (x axis)

Extrapolate the curve back to 3.0 minutes to establish the maximum temperature rise as shown in the example below:



Calculations:

1. Calculate the heat (q) released using the formula:

$$q = mc\Delta T$$

m = mass of solution (g)

c = specific heat capacity of water = $4.18 \text{ J g}^{-1} \text{ }^\circ\text{C}^{-1}$

ΔT = change in temperature ($^\circ\text{C}$)

2. Calculate the enthalpy change for the reaction between one mole of Zn and $\text{CuSO}_{4(\text{aq})}$

3. The literature value for the reaction is -217 kJ mol^{-1}

Calculate the percentage error using the following equation:

$$\% \text{ error} = \frac{\text{experimental value} - \text{literature value}}{\text{literature value}} \times 100$$

Title: Calculating the enthalpy change of neutralisation

Background: The enthalpy change of neutralisation (ΔH_n) is the enthalpy change when an acid and base react together to form one mole of water under standard conditions (100 kPa and 298 K). In this experiment you will calculate the enthalpy change of neutralisation between a strong acid hydrochloric acid (HCl) and strong base sodium hydroxide (NaOH).

Procedure:

1. Using a measuring cylinder, measure 30 cm³ of 1.00 moldm⁻³ HCl (hydrochloric acid) and pour into an insulated cup. Record the mass of the solution.
2. Insert a temperature probe into the cup and wait for the temperature to stabilise. Record the initial temperature in a table.
3. Using a measuring cylinder, measure 40 cm³ of 1.00 moldm⁻³ NaOH (sodium hydroxide). Measure the mass and temperature of the solution and record in a table.
4. Press collect on the data logger and the pour the 25cm³ of 1.00 moldm⁻³ NaOH into the insulated cup containing 25cm³ of 1.00 moldm⁻³ HCl
5. Collect the data for 180 seconds and record the maximum temperature reached.

Results:

Reactant	mass ($\pm 0.01\text{g}$)	initial temperature ($\pm 0.2\text{ }^\circ\text{C}$)	maximum temperature ($\pm 0.2\text{ }^\circ\text{C}$)	change in temperature ($\pm 0.4\text{ }^\circ\text{C}$)
HCl				
NaOH				

Calculations:




1. Determine the limiting reactant in the reaction
2. Calculate the heat released using $q = mc\Delta T$
3. Divide the heat released by the amount in mol of the limiting reactant.

Enthalpy change of neutralisation = _____ kJ mol⁻¹

Title: 6.1 Iodine clock reaction

To investigate the reaction of iodate(V) ions with sulfate(IV) ions in acidic solution and to determine the order of the reaction with respect to hydrogen ions.

Background: The iodate(V) ion is an oxidising agent and reacts with sulfate(IV) ions in acidic solution to produce iodine in solution, according to the following equation:



Iodine is only liberated when acid is added. It is possible to determine the effect of acid on the initial rate of reaction by timing how long it takes for the iodine to be produced. This is indicated by the formation of a blue colour with starch solution under different conditions of acid concentration, whilst keeping all other concentrations constant.

Apparatus:

300 cm³ beaker

Two 100 cm³ measuring cylinder

One 10 cm³ measuring cylinder

Potassium iodate(V) solution (200 cm³ 3.50 g dm⁻³)

Sodium sulfate(IV) solution (200 cm³ 5.00 g dm⁻³)

0.100 mol dm⁻³ sulfuric acid (500 cm³)

Starch solution (freshly prepared)

Stopwatch

Glass stirring rod

Procedure:

1. Measure 15 cm³ of the potassium iodate(V) solution into a 100 cm³ measuring cylinder and add 85 cm³ of deionised water to dilute it to 100 cm³.
2. Transfer this solution to a 400 cm³ beaker and add 5 cm³ of starch solution using a 10 cm³ measuring cylinder.
3. In a second 100 cm³ measuring cylinder, take 15 cm³ of sodium sulfate(IV) solution and add 85 cm³ of 0.100 mol dm⁻³ sulfuric acid.
4. Pour this solution into the 400 cm³ beaker containing the potassium iodate(V) solution and start the stop watch simultaneously. Stir the contents of the beaker with the glass rod.
5. Record in a Table the time taken (*t* in seconds) for the blue colour to appear.
6. Repeat the experiment using the same volumes of potassium iodate(V) solution, starch solution and sodium sulfate(IV) solution, but vary the concentration of sulfuric acid by changing the volume of sulfuric acid and making the total volume of the acidified sodium sulfate(IV) solution up to 100 cm³ by the addition of deionised water. The necessary volumes are shown in the table below with the initial experiment recorded on the right-hand side of the Table.

Experiment	1	2	3	4	5	Initial
Acid volume (cm ³)	25	35	45	55	70	85
Water volume (cm ³)	60	50	40	30	15	0
Na ₂ S ₂ O ₃ volume (cm ³)	15	15	15	15	15	15

Results:

Experiment	Concentration of H ₂ SO ₄ (mol dm ⁻³)	Time taken for blue colour to appear (s)
Initial		
1		
2		
3		
4		
5		

Topic 8.3 - Students should be familiar with the use of a pH meter and universal indicator.

Title: Reactions of acids and bases

Experiment 1 – Comparing the reaction of magnesium with a weak and strong acid

- Cut four 1 cm pieces of Mg ribbon
- Fill two test tubes with approximately 10cm³ of 1.0 mol dm⁻³ CH₃COOH and 1.0 mol dm⁻³ HCl
- Add two pieces of Mg ribbon to each test tube.
- Place a boiling tube over the top of the test tube to collect the gas produced.

Reaction	Observations
CH ₃ COOH + Mg	
HCl + Mg	

Experiment 2 – Comparing conductivity of a weak and strong acid and base

- Fill four 100 cm³ beakers with approximately 50 cm³ of 1.0 mol dm⁻³ CH₃COOH , 1.0 mol dm⁻³ HCl, 50 cm³ 1.0 moldm⁻³ NH₃ and 50 cm³ 1.0 mol dm⁻³ NaOH
- Set up an electrical circuit using graphite rods as electrodes.
- Measure the current that passes through each solution.

Electrolyte	Observations
1.0 mol dm ⁻³ CH ₃ COOH	
1.0 mol dm ⁻³ HCl	
1.0 mol dm ⁻³ NH ₃	
1.0 mol dm ⁻³ NaOH	

Experiment 3 – Comparing the pH of strong and weak acids and bases

- Use the same solutions as the previous experiment.
- Use a pH probe to test the pH of each of the solutions.
- Use pH paper to test the pH of each of the solutions.

Solution	pH with probe	pH with pH paper
1.0 mol dm ⁻³ CH ₃ COOH		
1.0 mol dm ⁻³ HCl		
1.0 mol dm ⁻³ NH ₃		
1.0 mo ldm ⁻³ NaOH		

Analysis:

Experiment 1

1) Write balanced chemical equations for the reactions of Mg with CH_3COOH and HCl .

2) Why does the strong acid react faster than the weak acid with Mg? (think about the factors that affect the rate of a reaction and the rate of formation of hydrogen gas).

3) Would you expect the same results if a weak and strong base were used? Explain your answer.

Experiment 2

4) What in the solutions allows the current to flow through?

5) Explain the difference in the conductivity of the weak and strong acids and bases.

Experiment 3

6) Explain the difference in pH between the strong and weak acids and bases.

7) Explain the difference between the strength and concentration of acids and bases.

19.1 - Perform lab experiments which could include single replacement reactions in aqueous solutions.

















Title: Displacement reactions

Background: Some metals are more reactive than others. The activity series lists metals in order of their reactivity (or strength as reducing agents). More reactive metals are at the top (stronger reducing agents) and less reactive metals are at the bottom (weaker reducing agents). In a displacement reaction, a more reactive metal displaces the ions of a less reactive metal from solution. In this experiment, you will react four metals with the salts of their ions.

The metals are zinc (Zn), lead (Pb), copper (Cu) and aluminium (Al). The metal salts are zinc sulfate, copper sulfate, lead sulfate and magnesium sulfate (all with a concentration of $0.100 \text{ mol dm}^{-3}$).

Procedure:

1. Add a few drops of each metal salt to the spotting tile as shown below.
2. Add a piece of each metal to the metal salts in the spotting tile.
3. Record your observations in the table below.

	Zn	Mg	Cu	Pb
ZnSO ₄				
MgSO ₄				
CuSO ₄				
PbSO ₄				

Results:

	Zn	Mg	Cu	Pb
ZnSO₄				
MgSO₄				
CuSO₄				
PbSO₄				

Analysis:

- 1) Arrange the metals in order of reactivity (from most reactive to least reactive).
Explain your reasoning.

- 2) Arrange the metals in order of the strength as reducing agents (from strongest to weakest).

- 3) Write net-ionic equations for the reactions that occurred.

9.2 Voltaic cells - Performance of laboratory experiments involving a typical voltaic cell using two metal/metal-ion half-cells.

Title: Construct a voltaic cell

Background: A voltaic cell is prepared from two half-cells connected together by a salt-bridge. A spontaneous reaction produces an electrical potential (voltage). A voltaic cell is constructed by using two metal electrodes in solutions of their respective salts (the electrolyte component of the cell) with known molar concentrations. In this experiment, you will use a voltmeter to measure the potential of a voltaic cell made with copper and zinc half-cells.

Procedure:

1. Use two 50 cm³ beakers to make the half cells.
2. Measure 20 cm³ of 0.10 mol dm⁻³ Cu(NO₃)₂ and 20 cm³ of 0.10 mol dm⁻³ Zn(NO₃)₂ in two 50 cm³ beakers.
3. Polish one Cu and one Zn metal strip for the electrodes.
4. Place the Cu strip in the beaker of Cu(NO₃)₂ solution and place the Zn strip in the beaker of Zn(NO₃)₂ solution. These are the half cells of your Cu-Zn voltaic cell.
5. Make a salt bridge by soaking a strip of filter paper in a beaker than contains a small amount of 1.00 mol dm⁻³ KNO₃ solution.
6. Connect the Cu and Pb half cells with the salt bridge.
7. Connect the wires to the electrodes and use a voltmeter to measure the potential of the Cu-Zn voltaic cell.

Analysis:

- 1) Draw a diagram of the voltaic cell below. Show the direction of electron flow.

- 2) Write half-equations that occur at the anode and cathode.